



Press Release

Optimum light and charge status indicator make for easier laryngoscopy

HEINE Standard F.O. L® LED Rechargeable Laryngoscope Handle / HEINE Standard F.O. L® Rechargeable Laryngoscope Handle

The new bottom cap, developed by HEINE for its laryngoscope handles, is exceptional. It now shows how much energy remains in the handle.

This happens fully automatically: the battery charge status is displayed on the handle for 10 seconds whenever it is removed from the charger.

The rechargeable laryngoscope handle completely recharges in just two hours in the NT300 – and features overcharge protection and no memory effect.

Compatible in every way: HEINE's development engineers always keep in mind that such an investment has to pay off in everyday use. As a result, it is possible to upgrade existing laryngoscope handles easily. Simply switch to a Li-ion L rechargeable battery and to the Standard F.O. L bottom cap.

The physician has the choice of two handles for the laryngoscope: either with a battery insert with the proven XHL Xenon Halogen bulb or with innovative LED illumination. Optimum lighting and comfortable working conditions are ensured no matter which he chooses.