

Combined Diagnostic Sets

COMBINED SETS

The HEINE Combined Diagnostic Sets contain three different instruments or more and include a range of high quality HEINE products.






HEINE Combined Diagnostic Sets




Set complete: **BETA 200 Ophthalmoscope, BETA 100 Diagnostic Otoscope**
 1 set of reusable specula
 10 disposable specula
 straight laryngeal mirror
 one spare bulb each
 hard case

Diagnostic Sets	3.5V (NiMH)	3.5V LI ION
BETA R rechargeable handle	A-231.20.376	A-231.27.376
BETA NT rechargeable handle and NT 300 charger	A-231.20.420	
 BETA L rechargeable handle and NT 300 charger		A-231.29.420

*Special Order - 4 to 6 weeks delivery time.




Set complete: **BETA 200 Ophthalmoscope, BETA 100 Diagnostic Otoscope**
 1 set of reusable specula
 10 disposable specula
 straight laryngeal mirror with accessories
 spreadable nasal speculum
 one spare bulb each
 zipper case

Diagnostic Sets	3.5V (NiMH)	3.5V LI ION
BETA R rechargeable handle	A-232.20.376	A-232.27.376
BETA NT rechargeable handle and NT 300 charger	A-232.20.420	
 BETA L rechargeable handle and NT 300 charger		A-232.29.420

*Special Order - 4 to 6 weeks delivery time.

As A-232 above, but with curved (instead of straight) laryngeal mirror with:



Diagnostic Sets	3.5V (NiMH)	3.5V LI ION
BETA R rechargeable handle	A-233.20.376	A-233.27.376
BETA NT rechargeable handle and NT 300 charger	A-233.20.420	
 BETA L rechargeable handle and NT 300 charger		A-233.29.420

*Special Order - 4 to 6 weeks delivery time.

: Rechargeable handle with charge status display in the bottom insert.

HEINE Combined Diagnostic Sets



Set complete: **BETA 200 Ophthalmoscope, BETA 200 F.O. Otoscope**
 1 set of reusable tips
 straight laryngeal mirror
 nasal illuminator with speculum
 Finoff transilluminator
 one spare bulb each
 hard case

Diagnostic Sets	3.5V (NiMH)	3.5V
BETA R rechargeable handle	A-157.20.376	A-157.27.376
BETA NT rechargeable handle and NT 300 charger	A-157.20.420	
BETA L rechargeable handle and NT 300 charger		A-157.29.420

*Special Order - 4 to 6 weeks delivery time.

: Rechargeable handle with charge status display in the bottom insert.

: To order your kit with the TL "Twist Lock" connection style, simply add the suffix "TL" to the end of the item number. I.E. B-141.20.376 TL for the set with Twist Lock.

HEINE mini 3000® Combined Diagnostic Sets



All handles are fitted with batteries (IEC LR 6 AA size, or optional rechargeable battery with mini NT charger*). Please order spare bulbs separately. The mini 3000 instruments in these sets can be supplied in black or blue. Please specify the color you prefer when ordering. If not otherwise specified, we will supply black instruments.

Set D-854 complete with **mini 3000 Ophthalmoscope, mini 3000 Otoscope, mini 3000 Combilamp**
 mini 3000 Laryngeal mirror
 1 set = 4 reusable tips
 hard case

mini 3000 Diagnostic Sets	2.5V
with mini 3000 battery handle with batteries	D-854.10.021
as D-854 but with nose speculum, 10mm dia.	D-853.10.021

*Special Order - 4 to 6 weeks delivery time.



To upgrade to a rechargeable instrument, with mini NT charger, please see page 146.

120V is the standard mains voltage. If a different voltage is required, e. g. 230V, please specify when ordering.



HEINE EN 100® Diagnostic Center



[01]



[02]

- ⚡ Soft Start Technology
- 🔒 Lockable handles
- ⚡ Automatic on/off switch
- ➕ Optional Extension Unit

An attractive modern design to keep your instruments in the right place, instantly ready for use.

With the optional Extension Unit as original equipment or for retrofitting to the EN 100 you have a system which lets you use three different instruments simultaneously. The EN 100 can also be combined with the AllSpec Tip Dispenser and the GAMMA XXL LF-W Sphygmomanometer to provide a complete diagnostic center to satisfy the most demanding user.

How to select your personal Diagnostic Center:

Order the EN 100 Wall Transformer and Extension Unit or additional system components such as the AllSpec Tip Dispenser or GAMMA XXL LF-W and your choice of 3.5V instruments individually from the HEINE range.

EN100 Diagnostic Center	AV	TL
EN 100* (Basic Unit with two handles) [01]	X-095.12.110	X-095.12.110 TL
EN 100-1* (Single Unit with one handle) [02]	X-095.12.109	X-095.12.109 TL
Extension Unit*	X-095.12.135	X-095.12.135 TL
AllSpec Tip dispenser [03]		B-000.11.148
GAMMA XXL LF-W Sphygmomanometer [04]		M-000.09.323

* without instruments



[05]

[04]

[01]


[03]

EN 100 Sets	AV	TL
with Wall Transformer EN 100 [01], BETA 200 F.O. Otoscope, BETA 200 Ophthalmoscope	A-095.12.976	A-095.12.976 TL
with Wall Transformer EN 100 [01], K180 F.O. Otoscope, K180 Ophthalmoscope	A-095.12.973	A-095.12.973 TL
with Wall Transformer EN 100 [01], K 180 F.O. Otoscope, K 180 Ophthalmoscope	A-095.12.974	A-095.12.974 TL
with Wall Transformer EN 100 [01], K 180 F.O. Otoscope, K 180 Ophthalmoscope, GAMMA XXL LF-W Sphygmomanometer [04], Tip Dispenser [03], on Wall Board [05].	A-095.12.700	A-095.12.700 TL
with Wall Transformer EN 100 [01], BETA 200 F.O. Otoscope, BETA 200 Ophthalmoscope, GAMMA XXL LF-W Sphygmomanometer [04], Tip Dispenser [03], on Wall Board [05].	A-095.12.701	A-095.12.701 TL
Diagnostic Center Wallboard [05]	Y-096.50.501	Y-096.50.501

A company makes history and gives a name to outstanding product quality: HEINE.

Founded in 1946, family owned, rich in tradition and today international market leader in a complete range of diagnostic instruments. HEINE products meet all international standards (ISO/CE) and represent the leading edge in precision and ergonomic design. Ongoing dedication to research and development has created a solid and diverse technology base for continued quality and product leadership. Over 500 employees worldwide contribute to this success.

www.heine.com



■ HEINE SUBSIDIARY COMPANIES
■ HEINE REPRESENTATIVES

NORTH AMERICA

HEINE USA LTD.
HEINE INSTRUMENTS CANADA LTD
10 Innovation Way
Dover, NH 03820
Tel. (603) 742-7103
Fax (603) 742-7217
Toll Free (800) 367-4872
E-Mail: service@heine-na.com

GERMANY

HEINE Optotechnik GmbH & Co. KG
82211 Herrsching
Kientalstr. 7
Tel. +49 (0) 81 52-38 0
Fax +49 (0) 81 52-3 82 02
E-Mail: info@heine.com

AUSTRALIA

HEINE AUSTRALIA PTY. LTD.
Unit 9, 98 Old Pittwater Road
P.O. Box 717 Brookvale 2100
New South Wales
Tel. +61 (0) 2-99 38 95 00
Fax +61 (0) 2-99 39 23 05
E-Mail: info@heine.com.au

SWITZERLAND

HEINE (Schweiz) AG
Tobeläckerstr. 9
CH-8212 Neuhausen am Rheinfall
Tel. +41 (0) 52-6 72 22 66
Fax +41 (0) 52-6 72 63 77
E-Mail: heine-schweiz@schaffhausen.ch

Your distributor

11/11. A-000.00.100 e

